

## Sterile Safety Syringe For Single Use



- Needle bevel oriented to the safety cover to facilitate low-angle injections

- Three point patent lock mechanism, once activated, the needle remains securely locked

- Easy to use, allow one-handed activation via two ways (thumb or surface)


- Audible and tactile click confirms safety activation

High-transparency syringe barrel

Latex-free rubber



Luer lock and Luer slip



## Syringe For Single Use With Safety Needle

Size	1ML,3ML,5ML,10ML,20ML,30ML,50ML and includes 30G, 29G, 28G, 27G, 26G, 25G, 24G, 23G, 22G, 21G, 20G, 19G, 18G.	Luer lock and Luer slip
Material	Barrel	Medical grade PP
	Plunger	Medical grade PP
	Rubber piston	Medical grade Rubber
	Needle Hub	Medical grade PP
	Needle protector	Medical grade PP
	Hub protector	Medical grade PP
	Cannula	SUS304
	Lubricant	Silicon oil
Intended use	The sterile hypodermic syringes for single use with safety needle is used by medical organizations to inject the liquid medicine into muscle of patients in the clinic process.	
Latex free	Yes	
DEHP/PVD/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	EN ISO 7886-1:2017	
Certification	CE, ISO13485, FDA	