

BERPU®

Blood Collection Set

- Wide, flexible wings for better control.
- Different color codes for easy identification.

- Diverse combinations available.



Sharp & smooth cannula edges
make penetration painless




Compatible with all
available holders.

CE 0123

STERILE EO

Technical Data Sheet



1. Rubber sleeve 2. Non-patient end needle 3. Needle hub 4. Tubing
5. Wing hub 6. Patient end needle 7. Needle protector

Blood Collection Set

Gauge	18G-27G	
Material	Wing Hub	Medical grade PVC
	Needle Protector	Medical grade PP
	Tubing	Medical grade PVC
	Cannula	SUS304
	Luer adaptor hub	Medical grade ABS
	Rubber sleeve	Medical grade rubber
Intended use	The blood collection set is applicable for blood collection from patients.	
Latex free	Customized	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO7864:2016	
Certification	CE, ISO13485	