


BERPU®

Insulin Syringes For Single Use



High transparency barrel, smooth gliding

Low dead space
-full utilisation of insulin,
minimum waste


Special silicon coating for advanced injection comfort.

Ultra fine needle, large insulin flow

CE 0123

STERILE	EO
---------	----

• Latex free rubber



Insulin Syringes

Size	0.3ml, 0.5ml, 1ml	
Material	Plunger	Medical grade PP
	Hub	Medical grade PP
	Hub protector	Medical grade PP
	Needle protector	Medical grade PP
	Rubber piston	Latex-free rubber
	Cannula	SUS304
	Lubricant	Silicon oil
Type/scale	U40 and U100, 0.3ml, 0.5ml, 1ml	
Intended use	The Insulin Syringe for single use is used for subcutaneous injection of insulin.	
Latex free	Yes	
DEHP/PVD/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO 8537:2016	
Certification	CE, ISO13485, FDA	

Technical Data Sheet



Unit Packaging: Both polybag and blister packaging are available

Packaging information for Insulin Syringe

Code	Insulin Unit	Volume	Needle size	Unit packaging Blister paper	Unit packing Polybag
I140052913	U-40	0.5ml	29GX1/2"	1pc/blister paper 100pcs/box 20 boxes/case, Csse dimension: 450x310x435mm	10pcs/polybag 100pcs/box 5 boxes/case, Csse dimension: 430x175x150mm
I140053008		0.5ml	30GX5/16"		
I140102913		1ml	29GX1/2"		
I140103013		1ml	30GX1/2"		
I140103108		1ml	30GX5/16"		
I110032913		U-100	0.3ml		
I110052913	0.5ml		29GX1/2"		
I110053013	0.5ml		30GX1/2"		
I110102913	1ml		29GX1/2"		
I110103008	1ml		30GX5/16"		
I110103108	1ml		30GX5/16"		

BERPU®

Tel: 86 577 85501388

E-mail: sales@berpu.com

Website: www.berpu.com

Address: NO.14 Xingji Road, Yongxing Street, Longwan District,
325000 Wenzhou Zhejiang Province, P.R. China