

BERPU[®]

Sterile Safety Pen Needles

For Single Use



Simple and
Beautiful appearance

- Ultra-fine needle
-special cut for gentle injection
- Large barrel
-easy handing and attachment
- Compatible with all major pens in the market.



Special silicon coating for
advanced injection comfort.

Patented automatic safety mechanism
-prevent needlestick

Before use

After use



CE 0123



STERILE EO



Safety Pen Needles

| | | |
|----------------------|--|--|
| Size | 29G to 31G, needle length from 4mm to 12mm | |
| Material | Needle protector | PP |
| | Needle base | PP |
| | Barrel | PP |
| | Spring | SUS304 |
| | Cannula | SUS304 |
| | Lubricant | Silicon oil |
| | Intended use | The Safety Pen Needles for single use are used with insulin pen for subcutaneous injection of insulin. |
| Latex free | Yes | |
| DEHP/PVD/BPA free | Yes | |
| Sterilization method | EO | |
| Shelf life | 5 years | |
| Standard | ISO 11608-2, EN ISO23908:2011 | |
| Certification | CE, ISO13485, FDA | |



Packaging:

Packaging information for Safety Pen Needles

| Code | Size | BOX QTY | CASE QTY | CASE DIMENSIONS |
|--------|----------|---------|----------|-----------------|
| I33104 | 31GX4MM | 100/box | 24/case | 380x325x318 |
| I33106 | 31GX6MM | 100/box | 24/case | 380x325x318 |
| I33108 | 31GX8MM | 100/box | 24/case | 380x325x318 |
| I33005 | 30GX5MM | 100/box | 24/case | 380x325x318 |
| I33006 | 30GX6MM | 100/box | 24/case | 380x325x318 |
| I33008 | 30GX8MM | 100/box | 24/case | 380x325x318 |
| I33010 | 30GX10MM | 100/box | 24/case | 380x325x318 |
| I32912 | 29GX12MM | 100/box | 24/case | 380x325x318 |

BERPU®

Tel: 86 577 85501388

E-mail: sales@berpu.com

Website: www.berpu.com

Address: NO.14 Xingji Road, Yongxing Street, Longwan District,
325000 Wenzhou Zhejiang Province, P.R. China