

BERPU®

Blood Collection Set

- Wide, flexible wings for better control.
- Different color codes for easy identification.

- Diverse combinations available.



Sharp & smooth cannula edges make penetration painless




Compatible with all available holders.

CE 0123

STERILE EO

Technical Data Sheet



1. Rubber sleeve 2. Non-patient end needle 3. Needle hub 4. Tubing
5. Wing hub 6. Patient end needle 7. Needle protector

Blood Collection Set

Gauge	18G-27G	
Material	Wing Hub	Medical grade PVC
	Needle Protector	Medical grade PP
	Tubing	Medical grade PVC
	Cannula	SUS304
	Luer adaptor hub	Medical grade ABS
	Rubber sleeve	Medical grade rubber
Intended use	The blood collection set is applicable for blood collection from patients.	
Latex free	Customized	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO7864:2016	
Certification	CE, ISO13485	

Technical Data Sheet



Packaging:

Packaging information for Blood collection set

Code	Inner packing	Size	BOX QTY	CASE QTY	DIMENSIONS
B92619	blister packing	26G×3/4"	100pcs/box	20/case	620*470mm*275mm
B92519	blister packing	25G×3/4"	100pcs/box	20/case	620*470mm*275mm
B92419	blister packing	24G×3/4"	100pcs/box	20/case	620*470mm*275mm
B92319	blister packing	23G×3/4"	100pcs/box	20/case	620*470mm*275mm
B92219	blister packing	22G×3/4"	100pcs/box	20/case	620*470mm*275mm
B92119	blister packing	21G×3/4"	100pcs/box	20/case	620*470mm*275mm
B92019	blister packing	20G×3/4"	100pcs/box	20/case	620*470mm*275mm
B92619	PE bag packing	26G×3/4"	100pcs/box	20/case	620*500mm*280mm
B92519	PE bag packing	25G×3/4"	100pcs/box	20/case	620*500mm*280mm
B92419	PE bag packing	24G×3/4"	100pcs/box	20/case	620*500mm*280mm
B92319	PE bag packing	23G×3/4"	100pcs/box	20/case	620*500mm*280mm
B92219	PE bag packing	22G×3/4"	100pcs/box	20/case	620*500mm*280mm
B92119	PE bag packing	21G×3/4"	100pcs/box	20/case	620*500mm*280mm
B92019	PE bag packing	20G×3/4"	100pcs/box	20/case	620*500mm*280mm

Tel: 86 577 85501388

E-mail: sales@berpu.com

Website: www.berpu.com

Address: NO.14 Xingji Road, Yongxing Street, Longwan District,
325000 Wenzhou Zhejiang Province, P.R. China

BERPU[®]