


BERPU®

Flashback Blood Collection Needle



Sharp & smooth cannula edges
-make penetration painless

Patented flashback design
-easier to operate, higher puncture
success rate




Thin-wall needle

-faster blood collection speed
and less hemolysis occurs

STERILE EO

CE 0123

Technical Data Sheet



1. Rubber sleeve 2. Bottom protector 3. Non-patient end needle
4. Non-patient end needle hub 5. Patient-end needle hub
6. Patient-end needle 7. Top protector

Flashback Blood Collection Needle

Size	20G to 22G, needle length from 25mm to 38mm	
Material	Needle Hub	Medical grade ABS
	Top Protector	Medical grade PP
	Bottom protector	Medical grade PP
	Cannula	SUS304
	Rubber sleeve	Medical grade rubber
	Lubricant	Silicone oil
Intended use	The blood collection needles for single use are applicable for blood collection from patients.	
Latex free	Customized	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO7864:2016	
Certification	CE, ISO13485	



Packaging:

Packaging information for Flashback Blood Collection Needle

Code	Size	BOX QTY	CASE QTY	DIMENSIONS
B22225	22GX1"	100/box	50/case	475x345x365
B22238	22GX1 1/2"	100/box	50/case	475x345x365
B22125	21GX 1"	100/box	50/case	475x345x365
B22138	21GX1 1/2"	100/box	50/case	475x345x365
B22025	20GX 1"	100/box	50/case	475x345x365
B22038	20GX1 1/2"	100/box	50/case	475x345x365

BERPU®

Tel: 86 577 85501388

E-mail: sales@berpu.com

Website: www.berpu.com

Address: NO.14 Xingji Road, Yongxing Street, Longwan District,
325000 Wenzhou Zhejiang Province, P.R. China