


BERPU[®]

Vacuum Blood Collecting Needle



Different color codes
for easy identification




Sharp & smooth cannula
edges make penetration
painless

Excellent self sealing rubber sleeve
ensures blood not flow out

CE 0123

STERILE EO



Vacuum Blood Collecting Needle

Size	16G to 23G, needle length from 25mm to 38mm	
Material	Needle Hub	Medical grade ABS
	Top Protector	Medical grade PP
	Bottom protector	Medical grade PP
	Cannula	SUS304
	Rubber sleeve	Medical grade rubber
	Lubricant	Silicone oil
	Intended use	The blood collection needles for single use are applicable for blood collection from patients.
Latex free	Customized	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO7864:2016	
Certification	CE, ISO13485, FDA	

Technical Data Sheet



Packaging:

Packaging information for Vacuum blood collecting needle

Code	Size	BOX QTY	CASE QTY	DIMENSIONS
B12325	23GX1"	100/box	50/case	415x365x295
B12225	22GX1"	100/box	50/case	415x365x295
B12238	22GX1 1/2"	100/box	50/case	415x365x295
B12125	21GX 1"	100/box	50/case	415x365x295
B12138	21GX1 1/2"	100/box	50/case	415x365x295
B12025	20GX 1"	100/box	50/case	415x365x295
B12038	20GX1 1/2"	100/box	50/case	415x365x295
B11838	18GX1 1/2"	100/box	50/case	415x365x295
B11638	16GX1 1/2"	100/box	50/case	415x365x295

BERPU®

Tel: 86 577 85501388

E-mail: sales@berpu.com

Website: www.berpu.com

Address: NO.14 Xingji Road, Yongxing Street, Longwan District,
325000 Wenzhou Zhejiang Province, P.R. China