



BERPU®

Sterile Safety Needles For Single Use

- Easy to use, allow one-handed activation via two ways (thumb or surface)
 - Audible and tactile click confirms safety activation
 - Needle bevel oriented to the safety cover to facilitate low-angle injections
 - Other sizes could be customized
- 




RW, TW, and ETW available



Match with various syringes with standard luer adapter

STERILE EO

CE 0123



Sterile Safety Needles For Single Use

Size	30G, 29G, 28G, 27G, 26G, 25G, 24G, 23G, 22G, 21G, 20G, 19G, 18G	
Material	Needle Hub	Medical grade PP
	Needle protector	Medical grade PP
	Hub protector	Medical grade PP
	Cannula	SUS304
	Lubricant	Silicon oil
Intended use	The Sterile safety hypodermic needles for single use are used for the injection of patients in the clinical process.	
Latex free	Yes	
DEHP/PVD/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	EN ISO 7864:2016	
Certification	CE, ISO13485, FDA	

Technical Data Sheet



Packaging:

Packaging information for Safety Needles

Code	Size	BOX QTY	CASE QTY	DIMENSIONS(mm)
N22716	27G*5/8"	100/box	20/case	510x415x235
N22620	26G*4/5"	100/box	20/case	510x415x235
N22520	25G*4/5"	100/box	20/case	510x415x235
N22425	24G*1"	100/box	20/case	510x415x235
N22325	23G*1"	100/box	20/case	510x415x235
N22332	23G*1 1/4"	100/box	20/case	510x415x235
N22232	22G*1 1/4"	100/box	20/case	510x415x235
N22138	21G*1 1/2"	100/box	20/case	510x415x235
N22038	20G*1 1/2"	100/box	20/case	510x415x235
N21938	19G*1 1/2"	100/box	20/case	510x415x235
N21838	18G*1 1/2"	100/box	20/case	510x415x235

BERPU®

Tel: 86 577 85501388

E-mail: sales@berpu.com

Website: www.berpu.com

Address: NO.14 Xingji Road, Yongxing Street, Longwan District,
325000 Wenzhou Zhejiang Province, P.R. China