

## Sterile Safety Syringe For Single Use



- Needle bevel oriented to the safety cover to facilitate low-angle injections

- Three point patent lock mechanism, once activated, the needle remains securely locked

- Easy to use, allow one-handed activation via two ways (thumb or surface)


- Audible and tactile click confirms safety activation

High-transparency syringe barrel

Latex-free rubber



Luer lock and Luer slip



## *Syringe For Single Use With Safety Needle*

|                      |   |                         |
|----------------------|---|-------------------------|
| Size                 | 1ML,3ML,5ML,10ML,20ML,30ML,50ML and includes 30G, 29G, 28G, 27G, 26G, 25G, 24G, 23G, 22G, 21G, 20G, 19G, 18G.   | Luer lock and Luer slip |
| Material             | Barrel  | Medical grade PP        |
|                      | Plunger   | Medical grade PP        |
|                      | Rubber piston   | Medical grade Rubber    |
|                      | Needle Hub  | Medical grade PP        |
|                      | Needle protector  | Medical grade PP        |
|                      | Hub protector   | Medical grade PP        |
|                      | Cannula   | SUS304                  |
|                      | Lubricant   | Silicon oil             |
| Intended use         | The sterile hypodermic syringes for single use with safety needle is used by medical organizations to inject the liquid medicine into muscle of patients in the clinic process. |                         |
| Latex free           | Yes   |                         |
| DEHP/PVD/BPA free    | Yes   |                         |
| Sterilization method | EO  |                         |
| Shelf life           | 5 years   |                         |
| Standard             | EN ISO 7886-1:2017  |                         |
| Certification        | CE, ISO13485, FDA   |                         |



## Packaging:

### Packaging information for Safety Syringes

| Code      | Size            | BOX QTY | CASE QTY | DIMENSIONS(mm) |
|-----------|-----------------|---------|----------|----------------|
| S6012325  | 1ML 23G*1"      | 100/box | 12/case  | 510x395x505    |
| S6032325  | 3ML 23G*1"      | 100/box | 12/case  | 510x395x475    |
| S6052230  | 5ML 22G*1 1/4"  | 100/box | 12/case  | 545x395x515    |
| S6102140  | 10ML 21G*1 1/2" | 100/box | 8/case   | 580x445x410    |
| S6201840  | 20ML 18G*1 1/2" | 100/box | 6/case   | 530x395x505    |
| S7012325  | 1ML 23G*1"      | 100/box | 12/case  | 510x395x505    |
| S7032325  | 3ML 23G*1"      | 100/box | 12/case  | 510x395x475    |
| S76052230 | 5ML 22G*1 1/4"  | 100/box | 12/case  | 545x395x515    |
| S7102140  | 10ML 21G*1 1/2" | 100/box | 8/case   | 580x445x410    |
| S7201840  | 20ML 18G*1 1/2" | 100/box | 6/case   | 530x395x505    |

# BERPU®

Tel: 86 577 85501388

E-mail: sales@berpu.com

Website: www.berpu.com

Address: NO.14 Xingji Road, Yongxing Street, Longwan District,  
325000 Wenzhou Zhejiang Province, P.R. China