


BERPU®

Dental Needles For Single Use

Smooth surface with light silicone coating
-minimal pain upon puncture

Universal plastic hub
-metric standard and British
standard threads, fit most syringes


Lancet bevel point at
cartridge end
-prevents blockage of
anesthetic



Triangle bevel mark
-for optimal injection position

STERILE EO

CE 0123



Dental Needles For Single Use

Size	30G, 27G, 25G	
Material	Needle tube	SUS304
	Needle hub	Medical grade PP
	Needle sheath	Medical grade PP
Intended use	The dental needles for single use is used with dental cartridge syringes for injection of dental local anaesthetics.	
Latex free	Yes	
DEHP/PVD/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO 7885-2010	
Certification	CE, ISO13485, FDA	

Technical Data Sheet



Packaging:

Packaging information for Blood Collection Tube

ITEM#	METRIC SIZE	BRITISH SIZE	BOX QTY	CASE QTY	DIMENSIONS(mm)
D13021	0.31x21mm	30Gx4/5"	100/box	50/case	395x315x302
D13025	0.31x25mm	30Gx1"	100/box	50/case	395x315x302
D12730	0.4x30mm	27Gx1 1/4"	100/box	50/case	395x315x302
D12735	0.4x35mm	27Gx1 2/5"	100/box	50/case	395x315x302
D12738	0.4x38mm	27Gx1 1/2"	100/box	50/case	395x315x302
D12530	0.5x30mm	25Gx1 1/2"	100/box	50/case	395x315x302

BERPU®

Tel: 86 577 85501388

E-mail: sales@berpu.com

Website: www.berpu.com

Address: NO.14 Xingji Road, Yongxing Street, Longwan District,
325000 Wenzhou Zhejiang Province, P.R. China