



# BERPU®

## Sterile Safety Needles For Single Use

- Easy to use, allow one-handed activation via two ways (thumb or surface)
  - Audible and tactile click confirms safety activation
  - Needle bevel oriented to the safety cover to facilitate low-angle injections
- 




RW, TW, and ETW  
available



Match with various syringes with  
standard luer adapter

STERILE EO

CE 0123



## *Sterile Safety Needles For Single Use*

Size	30G, 29G, 28G, 27G, 26G, 25G, 24G, 23G, 22G, 21G, 20G, 19G, 18G	
Material	Needle Hub	Medical grade PP
	Needle protector	Medical grade PP
	Hub protector	Medical grade PP
	Cannula	SUS304
	Lubricant	Silicon oil
Intended use	The Sterile safety hypodermic needles for single use are used for the injection of patients in the clinical process.	
Latex free	Yes	
DEHP/PVD/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	EN ISO 7864:2016	
Certification	CE, ISO13485, FDA	