

BERPU®



Sterile Hypodermic For Single Use

- Color marked polypropylene hubs
- Other sizes could be customized




- Self-made cannula for the production of hypodermic needle
- High flow
- Medical class stainless steel SUS304

- Regular wall, thin wall and extra thin wall are available
- Three bevel with sharp needle tip



CE 0123

STERILE EO



Hypodermic Needle For Single Use

Size	30G, 29G, 28G, 27G, 26G, 25G, 24G, 23G, 22G, 21G, 20G, 19G, 18G, 17G, 16G, 15G, 14G.	
Material	Needle Hub	Medical grade PP
	Needle protector	Medical grade PP
	Hub protector	Medical grade PP
	Cannula	SUS304
	Lubricant	Silicon oil
Intended use	The Sterile hypodermic needles for single use are used for the injection of patients in the clinical process.	
Latex free	Yes	
DEHP/PVD/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	EN ISO 7864:2016	
Certification	CE, ISO13485, FDA	